

Ellipse SuperLite Carbon Fibre Seat Walker

Lightweight walker suitable for indoor & outdoor use

- ✓ Foldable, lightweight Carbon Fibre Frame
- ✓ Height Adjustable Arms
- ✓ Added Cane Holder
- ✓ Includes Seat & Bag
- ✓ Curved Padded Backrest
- ✓ Maximum Weight Capacity – 120kg
- ✓ Aust. STD Tested 492484



Product Code	Colour	Size	Warranty (years)
9346376000635	Black	Small	1
9346376000611	Black	Medium	1
9346376000628	Black	Large	1

Specifications	Measurements
Length (cm)	76
Width (cm)	60
Maximum Height (cm)	105
Minimum Height (cm)	92
Folded Width (cm)	26
Seat Depth (cm)	22
Seat Width (cm)	45.4
Seat Height from the Ground (cm)	60
Front Castor Diameter (cm)	24.2/9.5"
Rear Castor Diameter (cm)	20.3/8"
Maximum Load (kg)	120
Net Weight (kg)	6.69