

Excel Superlite Manual Wheelchair

Lightweight manual wheelchair easy to transport

- ✓ Ultralite double braced aluminium frame
- ✓ Easily transportable - only 12.7kg
- ✓ Quick release rear wheels self-propel spoked (60cm) or transit (30cm)
- ✓ No inflation polyurethane "Soft Ride" tyres
- ✓ Foldable backrest for easy transportation and flip back padded armrest
- ✓ Swing away removable legrests
- ✓ Spill proof nylon fabric



Product Code	AWC133	AWC132
Colour	Blue	Blue
Seat Width (cm)	46	51
Seat Depth (cm)	43	43
Seat Height (cm)	50	50
Overall Width (cm)	64	69
Overall Length (cm)	107	107
Overall Height (cm)	95	95
Armrest Height (cm)	22	22
Backrest Height (cm)	45	45
Footrest Length (cm)	30 - 40	30 - 40
Folded Width (cm)	29	29
Rear Wheel Diameter (cm)	60	60
Front Castor Diameter (cm)	15	15
Max Load (kg)	110	110
Unit Weight (kg)	12.7	12.7
Warranty (years)	1	1

We care about
your comfort!

www.capeabilities.com.au

26 Strickland Street, Bunbury WA 6230
08 9717 6622
bunbury@capeabilities.com.au

Unit 1/64 West Street, Busselton WA 6280
08 9751 1446
sales@capeabilities.com.au