

Excel Superlite Transit Wheelchair

Lightweight transit wheelchair easy to transport




- ✓ Ultralite double braced aluminium frame
- ✓ Easily transportable - only 11.8kg
- ✓ Hand brakes (with park locking mechanism) - connected to drum hub brake system
- ✓ No inflation polyurethane "Soft Ride" tyres
- ✓ Foldable backrest for easy transportation and flip back padded armrest
- ✓ Swing away removable legrests
- ✓ Spill proof nylon fabric





Product Code	AWC136	AWC137
Colour	Blue	Blue
Seat Width (cm)	46	51
Seat Depth (cm)	43	43
Seat Height (cm)	50	50
Overall Width (cm)	64	69
Overall Length (cm)	92	92
Overall Height (cm)	95	95
Armrest Height (cm)	21	21
Backrest Height (cm)	45	45
Footrest Length (cm)	30 - 40	30 - 40
Folded Width (cm)	29	29
Rear Wheel Diameter (cm)	30	30
Front Castor Diameter (cm)	15	15
Max Load (kg)	110	110
Unit Weight (kg)	11.8	11.8
Warranty (years)	1	1

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au