


Funke Hyper Foam 2

The Funke Hyper Foam 2 mattress is a three-layer anti-decubitus therapy system.

Max weight capacity: 150kg






- Patient Comfort & Protection**

Ideal microclimate due to the hyper-elasticated foam mattress with longitudinal and diagonal intersections. Optimum pressure relief and minimum shear. PU-Tex cover is removable and washable up to 95 degrees. Also applicable for pain therapy and suitable for those up to category 4.

Product Code		Size
FUNKE-HF2-LS	Single	86 x 186 x 15cm
FUNKE-HF2-KS	King Single	105 x 196 x 15cm

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au