

ICare IC111

Adjustable bed including headboard

Maximum Weight Capacity (Long Single & King Single) 200kg

Maximum Weight Capacity (Double & Queen) 300kg






- The low-risk extra mobility option
- Offers head lift & knee break
- Fixed basic height (37cm)
- Lockable heavy duty castors
- External rechargeable battery back up
- End and corner retainer bars
- Hard wearing commercial grade fabric
- Strong breathable mattress platform
- Back lit wipeable hand piece
- No gap between mattresses on Queen & King beds
- Available in Long Single, King Single, Split Queen & Split King

Accessories Available :

- High side rail
- Low side rail
- Bed stick
- U-assist rail
- Over bed pole
- Full length side rails
- Headboard
- Footboard

**We care about
your comfort!**

www.capeabilities.com.au




 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au


 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au

Size	Long Single	King Single	Split Queen	Split King
Product Code	IC111LS	IC111KS	IC111Q	IC111K
Mattress Surface Size (mm)	2030 x 900	2030 x 1070	2030 x 750 (per side)	2030 x 900 (per side)
Base Size	2150 x 960	2150 x 1120	2150 x 780 (per side)	2150 x 900 (per side)
Hand Piece	Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable.		Individual hand control per side. Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable.	
Head Lift Angle	62 degrees			
Knee Break Angle	42 degrees			
Base Height	Fixed height of 37cm from floor to top of base.			
Actuators	1 x 6000 newton – 24 watt head adjust actuator. 1 x 6000 newton – 24 watt knee adjust actuator.			
Colour Choices	Onyx, Stone, Smoke, Silver, Vinyl Ebony (waterproof) & Vinyl Dove (waterproof).			
Warranty (against faulty workmanship and/or faulty materials)	5 years			

**We care about
your comfort!**

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au