

## ICare IC333

Function & durability. Features dual lift motors totalling 12,000 newtons lifting power, unlimited under-bed storage & 8 heavy duty castors

**Maximum Weight Capacity (Long Single & King Single) 200kg**

**Maximum Weight Capacity (Double & Queen) 300kg**



- Takes the hospital look out of the room
- Ultra low design reducing risk of falls
- Offers head lift & knee break
- Lockable heavy duty castors
- External rechargeable battery back up
- End and corner retainer bars
- Hard wearing commercial grade fabric
- Strong breathable mattress platform
- Back lit wipeable hand piece
- No gap between mattresses on Queen & King beds
- Available in Long Single, King Single, Long Double, One Piece Queen, Split Queen & Split King

### Accessories Available :

- High side rail
- Low side rail
- Bed stick
- U-assist rail
- Over bed pole
- Full length side rails
- Headboard
- Footboard

**We care about  
your comfort!**

[www.capeabilities.com.au](http://www.capeabilities.com.au)

26 Strickland Street, Bunbury WA 6230  
08 9717 6622  
[bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

Unit 1/64 West Street, Busselton WA 6280  
08 9751 1446  
[sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)

Size	Long Single	King Single	Long Double	One Piece Queen	Split Queen	Split King
Product Code	IC333LS	IC333KS	IC333D	IC333Q	IC333SQ	IC333K
Mattress Surface Size (mm)	2030 x 900	2030 x 1070	2030 x 1350	2030 x 1500	2030 x 750 (per side)	2030 x 900 (per side)
Base Size (mm)	2150 x 960	2150 x 1120	2150 x 1400	2150 x 1580	2150 x 780 (per side)	2150 x 900 (per side)
Hand Piece	Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable.		One hand control. Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable		Individual hand control per side. Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable.	One hand control. Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging on side of bed. Wipeable
Head Lift Angle	62 degrees					
Knee Break Angle	42 degrees					
Max Height	66cm from floor to top of base					
Min Height	22cm from floor to top of base					
Actuators	1 x 6000 newton – 24 watt head adjust actuator. 1 x 6000 newton – 24 watt knee adjust actuator. 2 x 6000 newton – 24 watt height adjust actuator (2 x 8000 newton on Double & Queen sizes).					
Colour Choices	Onyx, Stone, Smoke, Silver, Vinyl Ebony (waterproof) & Vinyl Dove (waterproof).					
Warranty (against faulty work-manship and/or faulty materials)	10 years					