

## Icare Mattresses | Firm, Medium & Soft

The ActiveX mattress offers support and comfort. This style is designed without springs, suitable for any profiling and especially suited to high risk users.

**Max weight capacity: 200kg**



- **3 Feels – Firm, Medium & Soft**

The **Firm Mattress** is ideal to help with mobility (easier turning).

The **Medium Mattress** is the most popular feel and provides support and mobility.

The **Soft Mattress** is most suited to those who spend a long time in bed (best suited to aged or petite patients).

- **Pressure Redistribution**

ActiveX is a high grade of elastic type foam which is responsive to temperature around inflamed areas. Around a pressure point, the mattress will soften to then allow blood to circulate freely.

- **XP Support Foam Base**




The solid support foam base is an open cell material that breathes to allow air circulation.

- **Patient Comfort & Protection**

The cover is made with the Bluetec cooling fabric, woven underneath with a cooling yarn to regulate body temperature. The purple vant fabric also has high air permeability and is anti-bacterial. For these beds, there are medical cover available (waterproof & wipe proof + welded seams & waterfall flaps).

We care about  
your comfort!

[www.capeabilities.com.au](http://www.capeabilities.com.au)




 26 Strickland Street, Bunbury WA 6230  
 08 9717 6622  
 [bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

 Unit 1/64 West Street, Busselton WA 6280  
 08 9751 1446  
 [sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)

Style	Size	IC15 (Firm) Depth - 15cm	IC20 (Medium) Depth - 20cm	IC25 (Soft) Depth - 22.5cm
Single	190 x 90cm	-	IC2S	ICXLS
Long Single	203 x 90cm	IC15LS	IC2LS	ICXLLS
King Single	203 x 107cm	IC15KS	IC2KS	ICXLKS
Double	190 x 135cm	-	IC2D	ICXLD
Long Double	203 x 135cm	IC15LD	IC2LD	ICXLLD
Queen	203 x 152cm	IC15Q	IC2Q	ICXLQ

**We care about  
your comfort!**

[www.capeabilities.com.au](http://www.capeabilities.com.au)

 26 Strickland Street, Bunbury WA 6230  
 08 9717 6622  
 [bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

 Unit 1/64 West Street, Busselton WA 6280  
 08 9751 1446  
 [sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)