

## Ibis Megan LRC – Standard & Small

Tilt in Space mechanism – Highly recommended by therapists for those with hip and lower back issues



- ✓ Tilt In Space mechanism
- ✓ Fully adjustable 2 or 3 pillow contoured backrest
- ✓ Available in 2 sizes – Standard & Small
- ✓ Chair is on castors for easier movement for cleaning
- ✓ Removable antimacassars and arm covers (fabric chairs only)
- ✓ Adjustable neckrest cushion & a removable seat protector
- ✓ Battery independent power source

Sizes	Standard	Small
Code	<b>CHMEG</b>	<b>CHMEGS</b>
Seat Height (cm)	51	48.5
Seat Width (cm)	53	53
Seat Depth (cm)	50	50
Backrest Height (cm)	70	68.5
Overall Width (cm)	74	74
<b>Max Load (kg)</b>	<b>142</b>	<b>142</b>
Warranty (years)	2	2

We care about  
your comfort!

[www.capeabilities.com.au](http://www.capeabilities.com.au)

26 Strickland Street, Bunbury WA 6230  
08 9717 6622  
[bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

Unit 1/64 West Street, Busselton WA 6280  
08 9751 1446  
[sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)