

Ibis LRC – Petite & Narrow

Customisable lift recliner







- ✓ Wallsaver mechanism allows chair to be positioned close to the wall
- ✓ Fully adjustable 2 or 3 pillow contoured backrest
- ✓ Available in 2 sizes – Standard & Narrow
- ✓ Chair is on castors for easier movement for cleaning
- ✓ Removable antimacassars and arm covers (fabric chairs only)
- ✓ Adjustable neckrest cushion & a removable seat protector
- ✓ Battery independent power source

Sizes	Standard (Petite)	Narrow
Code	CHPET	CHPETW
Seat Height (cm)	47	47
Seat Width (cm)	56	51
Seat Depth (cm)	47	47
Backrest Height (cm)	67	67
Overall Width (cm)	75	75
Max Load (kg)	142	142
Warranty (years)	2	2

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au