

Ikon Move Portable Travel Scooter




Portable scooter – Available in Metallic Blue or Metallic Orange



- ✓ Pulls apart for transport and storage
- ✓ Comfortable mid back contoured seat
- ✓ Full suspension
- ✓ Padded loop handlebar
- ✓ Angle adjustable tiller
- ✓ Width adjustable flip up padded armrests
- ✓ Removable shopping basket
- ✓ LED headlight
- ✓ Horn
- ✓ Fingertip operated controls
- ✓ Forward/reverse drive settings
- ✓ Variable speed setting
- ✓ Anti-tip wheels at rear



**We care about
your comfort!**

www.capeabilities.com.au




 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au



 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au

Product Code	IK4SEB	IK4SEO
Colour	Blue	Yellow
Overall length (mm)	1080	
Overall Width (mm)	590	
Overall Height (mm)	880	
Seat Width (mm)	450	
Seat Depth (mm)	440	
Seat Heights (mm)	330, 335, 380	
Armrest Length (mm)	320	
Armrest Width (mm)	80	
Turning Circle (mm)	1400	
Max Speed (km/h)	6.4	
Maximum Range (km)	Up to 18	
Safe Climb Angle	6°	
Maximum Climb Angle	10°	
Product Weight (kg)	60	
Maximum User Weight (kg)	130	

**We care about
your comfort!**

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au