

Jazzy Passport Compact Power Chair

Lightweight & effortless transportation

With convenient storage and foldable feature the Jazzy Passport makes travelling a breeze!


- ✓ Folds in a few simple steps - fold up & go!
- ✓ Back seat pocket and under seat mesh bag for extra storage.
- ✓ Includes a mesh cup holder and 101.6cm lap belt.
- ✓ Large foot platform to accommodate the user's comfort.
- ✓ Speciality foam seat design for extreme comfort.
- ✓ Dual braking system – regenerative & electromechanical
- ✓ Available with a 10.4 AH or 18 AH lithium ion air safe travel battery (check with your airline before flying).



Under-seat storage



Back seat pocket



Mesh cup holder



Product Features	Jazzy Passport
Weight Capacity	113kg
Maximum Speed	Up to 5.8km/h (3.6 mph)
Ground Clearance	7.9cm (3" at motor)
Turning Radius	6.16cm (24.25")
Overall Length	101.6cm (40")
Overall Width	59.7cm (23.5")
Seat to Ground Height	48.3cm (19")
Seat to Foot Platform Height	40cm (15.74") (no adjustment)
Casters	20.3cm (8") solid
Drive Tyres	30.5cm (12") pneumatic
Range Per Charge	Up to 15km (9.6 miles)
Standard Seating	Grey Foam, 45 x 43cm
Total Weight Without Batteries	24kg (54lbs)
Battery Weight	4kg (8.8lbs)

**We care about
your comfort!**

www.capeabilities.com.au

📍 26 Strickland Street, Bunbury WA 6230
📞 08 9717 6622
✉ bunbury@capeabilities.com.au

📍 Unit 1/64 West Street, Busselton WA 6280
📞 08 9751 1446
✉ sales@capeabilities.com.au