

Karma Ergo Lite Deluxe Manual Wheelchair

Ultra-lightweight manual wheelchair easy to transport

Two size options to choose from

- ✓ Ultra lightweight aluminium folding frame with double cross brace for strength
- ✓ Easily transportable only 10.2-10.3kg
- ✓ Patented S-shaped ergonomic seating system which promotes pressure relief and helps prevent the user sliding forward
- ✓ Ergonomic arm rests and height adjustable footrests
- ✓ Fire retardant upholstery seat and folding back
- ✓ Optional fixed or quick release 20" self-propel wheel
- ✓ Attendant propelled option available



Product Code	2-KM2512F20-16	2-KM2512F20-18
Colour	Silver	Silver
Total Height (cm)	91	91
Total Length (cm)	93	93
Total Width (cm)	57	63
Seat Width (cm)	41	45
Seat Depth (cm)	42	42
Seat Height (cm)	48	48
Armrest Height (cm)	21.5	21.5
Backrest Height (cm)	42	42
Folded Width (cm)	31	31
Max Load (kg)	100	100
Net Weight (kg)	10.2	10.3
Warranty (years)	1	1