

Karma Ergo Lite Transit Wheelchair

Ultra-lightweight transit wheelchair easy to transport

Two size options to choose from

- ✓ Ultra lightweight aluminium folding frame with double cross brace for strength
- ✓ Easily transportable – only 8.1-8.3kg
- ✓ Patented S-shaped ergonomic seating system which promotes pressure relief and helps prevent the user sliding forward
- ✓ Ergonomic arm rests and height adjustable footrests
- ✓ Fire retardant upholstery seat and folding back
- ✓ Self propelled option available



Product Code	2-KM2501-16	2-KM2501-16
Colour	Silver	Silver
Total Height (cm)	91	91
Total Length (cm)	94	94
Total Width (cm)	57	62
Seat Width (cm)	41	46
Seat Depth (cm)	41	41
Seat Height (cm)	48	48
Armrest Height (cm)	21.5	21.5
Backrest Height (cm)	42	42
Folded Width (cm)	27.5	27.5
Max Load (kg)	100	100
Net Weight (kg)	8.1	8.3
Warranty (years)	1	1

We care about
your comfort!

www.capeabilities.com.au

26 Strickland Street, Bunbury WA 6230
08 9717 6622
bunbury@capeabilities.com.au

Unit 1/64 West Street, Busselton WA 6280
08 9751 1446
sales@capeabilities.com.au