

## Let's Go Out Rollator

Lightweight, smooth & steady

- ✓ Easy to fold and will stand by itself when folded.
- ✓ Elegant large front wheels make it easy to turn and cross pavement curbs. Rear wheel fenders protect clothing from dirt.
- ✓ Shock-absorbing suspension on the rear wheels make it run smoothly on all surfaces.
- ✓ Ergonomically shaped rubber hands.
- ✓ Effective brake and parking brake, which are easy to operate with no visible wires.
- ✓ Seat is made of black/beige artificial leather. Soft and comfortable to sit on, as well as durable and water resistant.



Product Code	Colour	Size	Warranty (years)
0120	Black	One Size	1
0150	Beige	One Size	1

Specifications	Measurements
Width (cm)	58
Length (cm)	65
Weight (kg)	6.4
Handle Height is Adjustable from (cm)	78 - 96
Seat Height (cm)	59.5
Seat Width (cm)	46
Front Wheel Diameter (cm)	24.5
Rear Wheel Diameter (cm)	20
Max. Bag Load (kg)	10
Maximum User Weight (kg)	130

We care about  
your comfort!

[www.capeabilities.com.au](http://www.capeabilities.com.au)

26 Strickland Street, Bunbury WA 6230  
08 9717 6622  
[bunbury@capeabilities.com.au](mailto:bunbury@capeabilities.com.au)

Unit 1/64 West Street, Busselton WA 6280  
08 9751 1446  
[sales@capeabilities.com.au](mailto:sales@capeabilities.com.au)