

Shoprider Little Ripper Scooter

Small portable scooter suitable for indoor and outdoor usage

Two-colour options to choose from




- ✓ Easy to pull apart for storage and transport
- ✓ Flat free tyres
- ✓ Comfortable mid back contoured seat
- ✓ Removable/rotating/height adjustable seat.
- ✓ Lift up & adjustable width armrests.
- ✓ Removable shopping basket.



Product Code	GK9BL	GK9RE
Colour	Blue	Burgundy Red
Size	One Size	One Size
Overall Length (cm)	105	105
Overall Width (cm)	56	56
Overall Height (cm)	88	88
Turning circle (cm)	115	115
Max Speed (km/h)	6	6
Max Travel Range (km)	20	20
Safe Climb Angle (°)	9°	9°
Max Climb Angle (°)	9°	9°
Motor Size	0.8Hp (2 pole motor)	0.8Hp (2 pole motor)
Batteries	2 x 21 Amp (supplied)	2 x 21 Amp (supplied)
Charger	5 Amp	5 Amp
Wheel Size (cm)	23 (front & rear)	23 (front & rear)
Max Load (kg)	136	136
Product Weight (kg)	57.6 with batteries	57.6 with batteries
Battery Weight (kg)	7.7 each	7.7 each
Warranty (years)	1	1

We care about
your comfort!

www.capeabilities.com.au

 26 Strickland Street, Bunbury WA 6230
 08 9717 6622
 bunbury@capeabilities.com.au

 Unit 1/64 West Street, Busselton WA 6280
 08 9751 1446
 sales@capeabilities.com.au