

U-Step Walker 2

Created for people with Parkinsons to increase independence

- ✓ Stable and secure solid frame.
- ✓ Patented U-shaped base. Its ultra-stable foundation braces you in every direction.
- ✓ Rear wheel is located in the centre to your gravity and prevent falls. Also helps travel through tight spaces, doorways etc.
- ✓ Easy to use and effective braking system. Lightly squeeze the brake release lever to roll. Stable stationary frame when transferring from sitting to standing.
- ✓ Easy to fold in 2 steps – no need to bend down.
- ✓ Optional laser and sound cue extras to normalise and guide your walking.



Specifications	Measurements
Width (cm)	58
Length (cm)	63.5
Weight (with seat and basket) (kg)	10
Seat Height (cm)	56
Seat Width (cm)	51
Turning circle (cm)	74
Folded Dimensions (cm)	107 x 58 x 25
Patient Height Range (cm)	150 x 187
Patient Height Range (with customisation) (cm)	150 x 193
Maximum User Weight (kg)	170