

### Cleaning and drying your accessories

Remove the DnaBean from your accessory before cleaning.

For leather accessories:

1. Wipe them clean with a non-abrasive, lint-free cloth\*.
2. Do not soak them in water. They are not showerproof or water resistant. Soaking in water may cause some colour release.
3. Allow them to air-dry thoroughly before you reinsert the DnaBean. Do not store them in direct sunlight, at high temperatures or in high humidity.

For all other accessories:

1. Wipe them clean with a non-abrasive, lint-free cloth. If necessary, lightly dampen the cloth with fresh water.
2. Dry them with a non-abrasive, lint-free cloth before you reinsert the DnaBean.

\* Leather usually takes on a patina after normal wear. The leather may also change colour due to exposure to chemical substances like perfume, insect repellent, oil, lotion or dyed materials like denim. The colour of an accessory may vary or fade over time.



DN-ENG-IFU-017 Rev 0  
9501019 10-2021