
























1 | Intended Use

The device is intended to express milk from lactating women in order to collect milk from their breasts. Annabella is a home use device, but can also be used in other places (or similar environments). Annabella is a single electric breast pump, and can be used only on one breast at a time. The pump does not treat or is meant to treat inverted or sunken nipples.

2 | Meaning of Symbols

SYMBOL	DESCRIPTION
	Manufacturer
	Date of manufacture
	Batch code
	Catalogue number
	Serial number
	Fragile, Handle with care
	Keep away from sunlight
	Keep dry
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation

SYMBOL	DESCRIPTION
	Caution
	Both direct and alternating current
	Read and follow the instructions for use
	CLASS II equipment
	"ON" / "OFF" (push-push)
	TYPE BF APPLIED PART
IP22	Protected from touch by fingers and objects greater than 12 millimeters and/or Protected from water spray less than 15 degrees from vertical.
	The packaging contains products intended to come in contact with food
	Indoor Use Only
	The package is capable of being recycled (including washable parts)
UDI	Unique device identifier / Global Unique Device Identification Database
	Country of Origin
Patent Number	Patent Number
	Do not throw in trash. Send to a recycling center (electrical components)
BPA Free	BPA Free
	Lithium Battery Caution