

Warning!

DO NOT sanitize the parts in the microwave!




Clean

10 Min.




Dishwasher Safe




Important notice: the silicone and plastic parts have an expiration date and should be maintained and replaced approximately **every 3 months**.

Which parts should I sanitize?

Below is a list of the parts that can and should be sanitized before your first use:



Pumping with Annabella should be pleasant. Feel any pain? We are here to help!




Annabella

Quick User Guide

Assembly


1 | Charge




2 | Assemble




3 | Place




4 | Connect




5 | Push




6 | Insert




7 | Hold and Click




8 | Hold and Turn



9 | Connect



Check Your Fit




Just Right

Too Big

Too Small

Align




Pumping Session

Unplug from power source before pumping!


1 | Power On




2 | Start




3 | Align



5 | Vacuum (+/-)



6 | Tongue Pace (+/-)



7 | Tongue Height

