

Item #:	Item Name:	Description:	Product Image
<b>CC01117-T (Tris)</b> <b>CC01117-P (PBS)</b> *choose one	<b>FluidPrep EasyElute Large Volume Concentration Kit</b> *Requires a vacuum pump for operation	1. EasyElute™ 30 kD UltraFilter 2. EasyElute™ Elution Buffer canister (Tris or PBS) 3. Vacuum Tube (12") 4. Sample Tube with fitting (36") 5. Silicone Port Caps (x2) 6. Canister Adapter 7. Original Plastic Port Covers (x2) <b>These are discarded before use)</b>	
Item #:	Item Name:	Description:	Product Image
<b>CC02000-T (Tris)</b> <b>CC02000-P (PBS)</b> *choose one	<b>FluidPrep EasyElute Large Volume Concentration Kit</b> <b>With AUTOCLAVED TUBING</b> *Requires a vacuum pump for operation	1. EasyElute™ 30 kD UltraFilter 2. EasyElute™ Elution Buffer canister (Tris or PBS) 3. Vacuum Tube (12") 4. Sample Tube with fitting (36") 5. Silicone Port Caps (x2) 6. Canister Adapter 7. Original Plastic Port Covers (x2) <b>These are discarded before use)</b>	
Item #:	Item Name:	Description:	Product Image
<b>CC01118-T (Tris)</b> <b>CC01118-P (PBS)</b> *choose one	<b>FluidPrep EasyElute Resupply Kit</b>	1. EasyElute™ 30 kD UltraFilter 2. EasyElute™ Elution Buffer canister (Tris or PBS)	
Item #:	Item Name:	Description:	Product Image
<b>HC08018-T (Tris)</b>	<b>FluidPrep EasyElute LVC High Volume Elution Buffer</b>	Single-use Elution Buffer Canister for the EasyElute LVC. Please select PBS or Tris  The EasyElute Canister Adaptor (CC08102) is required for use.  *Do not use with CP Select, may cause damage.	
Item #:	Item Name:	Description:	Product Image
<b>HC08018-P (PBS)</b>	<b>FluidPrep EasyElute LVC High Volume Elution Buffer</b>	Single-use Elution Buffer Canister for the EasyElute LVC. Please select PBS or Tris  The EasyElute Canister Adaptor (CC08102) is required for use.  *Do not use with CP Select, may cause damage.	
Item #:	Item Name:	Description:	Product Image
<b>CC08102</b>	<b>EasyElute Canister Adaptor</b>	The canister adaptor is a replacement part for the EasyElute LVC Kit.	