

Product datasheet



Nine Health Latexgloves, white

Nine Health Latex gloves are powder-free and chlorinated latex gloves that are suitable for general care and cleaning tasks as well as food handling. The gloves are slightly nubbed (textured) for a better grip. The gloves contain natural latex, which for some can cause allergic reactions, including anaphylactic shock.



Size	XS	S	M	L	XL
Item no.	100120	100121	100122	100123	100124
Pcs/box	100	100	100	100	80

Material	Latex, powderfree, chlorinated
CE classification	MDR 2017/745 class 1
	PPE (EU) 2016 / 425 CAT III
	EN ISO 21420
Tests & quality Standards	EN 455-1, EN 455-2, EN 455-3, EN 455-4
	EN ISO 374-1, EN374-2, EN 374-4, EN ISO 374-5, EN ISO 16523-1
	EN 1186 Food
	EN 388
	Food Commission Regulation (EU) No. 10/2011



Product datasheet




Nine Health Latexgloves, white

Dimensions						
Size	Weight, G.	Lenght min.	Width - mm	Thickness - Manchet	Thickness - Palm	Thickness - Finger
X-Small	5,8	240 mm	75 ± 5	min. 0.10 mm	min. 0.12 mm	min. 0.14 mm
Small	6,3	240 mm	85 ± 5	min. 0.10 mm	min. 0.12 mm	min. 0.14 mm
Medium	6,8	240 mm	95 ± 5	min. 0.10 mm	min. 0.12 mm	min. 0.14 mm
Large	7,3	240 mm	106 ± 5	min. 0.10 mm	min. 0.12 mm	min. 0.14 mm
X-Large	7,8	240 mm	116 ± 5	min. 0.10 mm	min. 0.12 mm	min. 0.14 mm

AQL	1,5
Breakforce min.	7N before aging, 6N after aging
Shelf life	3 years
Elongation	650% before aging, 500% after aging
Keep dark and by roomtemperature	

Easily recognizable colored boxes indicating sizes - Also across product categories



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com