

Product datasheet



NINE UNITED CARE

Nine Health Nitril examination gloves, BerryBlue, powderfree

Non-sterile, blue coloured, low dermatitis, ambidextrous, finger-textured, on-line chlorinated powder-free nitrile examination glove. Tested to protect against Fentanyl for at least 240 minutes. Suitable for people with sensitive skin / allergies.



Size	XS	S	M	L	XL
Item no.	100070	100071	100072	100073	100074

Material	Nitril, powderfree
CE classification	MDR (EU)2017/745 Class I
	PPER (EU) 2016/425 Cat III
	FCM (EU) 1935/2004
Tests & quality Standards	EN 455-1, EN 455-2, EN 455-4, EN 455-5
	EN 420
	EN ISO 374-1, EN374-2, EN 374-4, EN ISO 374-5, EN ISO 16523-1
	EN ISO 13485
	EN ISO 14971
	Resolution ResAP (2004) 4 Food contact Material - Rubber
	ASTM D6978-05 (Chemoterapy drugs)*
	ASTM D7558-09 (Accelator)*
	* Test report available on request



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com

Product datasheet



Nine Health Nitril examination gloves, BerryBlue, powderfree

Dimensions						
Size	Weight, G.	Lenght min.	Width - mm	Thickness - Manchet	Thickness - Palm	Thickness - Finger
X-Small	2,8	240 mm	75 ± 5	min. 0.04 mm	min. 0.05mm	min. 0.07 mm
Small	3	240 mm	85 ± 5	min. 0.04 mm	min. 0.05mm	min. 0.07 mm
Medium	3,2	240 mm	95 ± 5	min. 0.04 mm	min. 0.05mm	min. 0.07 mm
Large	3,4	240 mm	106 ± 5	min. 0.04 mm	min. 0.05mm	min. 0.07 mm
X-Large	3,6	240 mm	116 ± 5	min. 0.04 mm	min. 0.05mm	min. 0.07 mm

AQL	1,5
Breakforce min.	6N
Shelf life	3 years
Elongation	500% before aging, 400% after aging
Keep dark and by roomtemperature	

Easily recognizable colored boxes indicating sizes - Also across product categories



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com