

# Product datasheet



## Nine Health Nitril examination gloves, black

Non-sterile disposable nitrile gloves that can be used for tasks where you are in contact with blood and other body fluids, handling medications, cytostatics, food handling, cleaning tasks, care tasks, etc.



<b>Size</b>	XS	S	M	L	XL
<b>Item no.</b>	100050	100051	100052	100053	100054

<b>Material</b>	Nitril, powderfree
<b>CE classification</b>	MDR (EU)2017/745 Class I
	PPER (EU) 2016/425 Cat III
	FCM (EU) 1935/2004
<b>Tests &amp; quality Standards</b>	EN 455-1, EN 455-2, EN 455-4, EN 455-5
	EN 420
	EN ISO 374-1, EN374-2, EN 374-4, EN ISO 374-5, EN ISO 16523-1
	EN ISO 13485
	EN ISO 14971
	Resolution ResAP (2004) 4 Food contact Material - Rubber



# Product datasheet



## Nine Health Nitril examination gloves, black

Dimensions						
Size	Weight, G.	Lenght min.	Width - mm	Thickness - Manchet	Thickness - Palm	Thickness - Finger
X-Small	3.3 ± 0.3	240 mm	75 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Small	3.7 ± 0.3	240 mm	85 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Medium	4.2 ± 0.3	240 mm	95 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Large	4.6 ± 0.3	240 mm	105 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
X-Large	5.0 ± 0.3	240 mm	115 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm

AQL	1,5
Breakforce min.	6N
Shelf life	3 years
Elongation	500% before aging, 400% after aging
Keep dark and by roomtemperature	

Easily recognizable colored boxes indicating sizes - Also across product categories



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

[www.nineunitedcare.com](http://www.nineunitedcare.com)