

Product datasheet




Nine Health Nitril examination gloves, blue, 12 inch / 29 cm

Non-sterile disposable nitrile gloves with extra long cuff that can be used for tasks where you are in contact with blood and other body fluids, handling medications, cytostatics, food handling, cleaning tasks, care tasks, etc.



Size	XS	S	M	L	XL
Item no.	100165	100166	100167	100168	100169

Material	Nitril, powderfree
CE classification	MDR (EU)2017/745 Class I
	PPER (EU) 2016/425 Cat III
	FCM (EU) 1935/2004
Tests & quality Standards	EN 455-1, EN 455-2, EN 455-4, EN 455-5
	ISO 21420
	ISO 374-1, ISO 374-2, ISO 374-4, ISO ISO 374-5, EN 16523-1
	EN ISO 13485
	EN ISO 14971
	Resolution ResAP (2004) 4 Food contact Material - Rubber
	ASTM D6978-05 (Chemoterapy drugs)*
	ASTM D7558-09 (Accelerator)*
	* Test report available on request



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com

Product datasheet



Nine Health Nitril examination gloves, blue, 12 inch / 29 cm

Dimensions						
Size	Weight, G.	Lenght min.	Width - mm	Thickness - Cuff	Thickness - Palm	Thickness - Finger
X-Small	6.0 ± 0.3	290 mm	75 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Small	6.5 ± 0.3	290 mm	85 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Medium	7.0 ± 0.3	290 mm	95 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
Large	7.4 ± 0.3	290 mm	105 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm
X-Large	7.8 ± 0.3	290 mm	115 ± 5	min. 0.05 mm	min. 0.07 mm	min. 0.11 mm

AQL	1,5
Breakforce min.	6N
Shelf life	3 years
Elongation	500% before aging, 400% after aging
Keep dark and by roomtemperature	

Easily recognizable colored boxes indicating sizes - Also across product categories



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com