

Product datasheet



Nine Health Vinyl examination gloves, transparent

Non-sterile disposable vinyl gloves contains no latex protein,
Nonallergic.

Pinhole rate is low to avoid contact with
infectious substances. They are ambidextrous,
wear convenient, comfortable with manual
dexterity. Approved for use in food handling.



Size	XS	S	M	L	XL
Item no.	100160	100161	100162	100163	100164

Material	Vinyl, PVC without DEHP, powderfree, latexfree
CE classification	MDR (EU)2017/745 Class I
Tests & quality Standards	EN 455-1, EN 455-2, EN 455-4, EN 455-5 ISO 21420 ISO 374-1, ISO 374-2, ISO 374-4, ISO ISO 374-5, EN 16523-1 ISO 13485 EN 13130-4 Resolution ResAP (2004) 4 Food contact Material - Rubber Regulation (EC) No 1907/2006 concerning REACH Regulation (EU) No 2020/1245 + No 10/2011 regarding food contact



Product datasheet



Nine Health Vinyl examination gloves, transperant

Dimensions						
Size	Weight, G.	Lenght min.	Width - mm	Thickness - Cuff	Thickness - Palm	Thickness - Finger
X-Small	3.5 ± 0.3	240 mm	75 ± 5	min. 0.04 mm	min. 0.06 mm	min. 0.07 mm
Small	4,0 ± 0.3	240 mm	85 ± 5	min. 0.04 mm	min. 0.06 mm	min. 0.07 mm
Medium	4,5 ± 0.3	240 mm	95 ± 5	min. 0.04 mm	min. 0.06 mm	min. 0.07 mm
Large	5,0 ± 0.3	240 mm	105 ± 5	min. 0.04 mm	min. 0.06 mm	min. 0.07 mm
X-Large	5.5 ± 0.3	240 mm	115 ± 5	min. 0.04 mm	min. 0.06 mm	min. 0.07 mm

AQL	1,5
Breakforce min.	3,6N
Shelf life	3 years
Keep dark and by roomtemperature	

Easily recognizable colored boxes indicating sizes - Also across product categories



NINE UNITED CARE A/S

Havnen 1, 8700 Horsens

www.nineunitedcare.com