

Denoix Series

- Sizes 5 -10 made from 3/8" thick aluminum (7 – 10 also available in 1/2")
- Sizes 11 – 14 made from 1/2" thick aluminum



Reverse Front

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 5/8	117.5	12 3/16	322.2
6	4 ¾	120.6	4 7/8	123.8	13 3/8	340.1
7	5	127	5 1/8	130.1	14 1/16	358
8	5 ¼	133.3	5 3/8	136.5	14 7/8	375.9
9	5 ½	139.7	5 11/16	144.4	15 ½	393.7
10	5 ¾	146	5 15/16	150.8	16 3/16	411.7
11	6	152.4	6 3/16	157.2	17	431.8
12	6 ¼	158.7	6 7/16	163.5	17 5/8	447.7
13	6 ½	165.1	6 11/16	170	18 3/8	466.7
14	6 ¾	171.5	6 15/16	176.2	19	482.6



Suspensory Front

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 15/16	125.4	12 7/8	327.6
6	4 ¾	120.6	5 3/16	131.7	13 1/8	345.9
7	5	127	5 7/16	138.1	14 5/16	364.2
8	5 ¼	133.3	5 3/4	146	15 1/16	382.3
9	5 ½	139.7	6	152.4	15 3/4	401
10	5 ¾	146	6 1/4	158.7	16 1/2	419
11	6	152.4	6 9/16	166.7	17 1/4	438.1
12	6 ¼	158.7	6 7/8	174.6	17 15/16	455.6
13	6 ½	165.1	7 1/16	179.4	18 5/8	473
14	6 ¾	171.5	7 3/8	187.3	19 3/8	492.1



Suspensory Hind

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	5	127	12 7/8	327.6

6	4 ¾	120.6	5 ¼	133.3	13 5/8	345.7
7	5	127	5 ½	139.7	14 5/16	363.7
8	5 ¼	133.3	5 13/16	147.6	15 1/16	382
9	5 ½	139.7	6 1/16	154	15 3/4	400
10	5 ¾	146	6 3/8	162	16 7/16	418.3
11	6	152.4	6 5/8	168.3	17 7/16	443
12	6 ¼	158.7	6 7/8	174.6	18 1/8	460
13	6 ½	165.1	7 1/8	181	18 7/8	479.4
14	6 ¾	171.5	7 7/16	189	19 5/8	498.5

Suspensory Branch (Hind)  and Collateral Hind 

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 15/16	125.4	13 1/16	331.2
6	4 ¾	120.6	5 3/8	136.5	13 3/4	349.7
7	5	127	5 1/2	139.7	14 1/2	368
8	5 ¼	133.3	5 3/4	146	15 3/16	386.5
9	5 ½	139.7	6 1/16	154	15 15/16	404.8
10	5 ¾	146	6 5/16	160	16 11/16	423.4
11	6	152.4	6 9/16	166.7	17 3/8	441.3
12	6 ¼	158.7	6 7/8	174.6	18 1/8	460
13	6 ½	165.1	7 1/8	181	18 7/8	480
14	6 ¾	171.5	7 7/16	190	19 1/2	495

Collateral Front 

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 5/8	117.5	12 3/4	323.3
6	4 ¾	120.6	4 7/8	123.8	13 7/16	341.3
7	5	127	5 1/8	130	14 1/8	359.4
8	5 ¼	133.3	5 3/8	136.5	14 7/8	377.4
9	5 ½	139.7	5 5/8	142.8	15 9/16	395.2
10	5 ¾	146	5 7/8	149.2	16 1/4	453.1
11	6	152.4	6 1/8	155.6	17	431.8
12	6 ¼	158.7	6 7/16	163.5	17 11/16	449.2
13	6 ½	165.1	6 11/16	170	18 3/8	467

14	6 ¾	171.5	6 15/16	176.2	19 1/16	484.2
----	-----	-------	---------	-------	---------	-------



Onion & Onion Lateral Roller

Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 ½	114.3	12 1/2	318.2
6	4 ¾	120.6	4 13/16	122.2	13 3/16	335.7
7	5	127	5 1/16	128.5	13 7/8	353.5
8	5 ¼	133.3	5 5/16	135	14 5/8	371.1
9	5 ½	139.7	5 5/8	142.8	15 5/16	388.8
10	5 ¾	146	5 7/8	149.2	16	406.6
11	6	152.4	6 1/8	155.6	16 11/16	423.8
12	6 ¼	158.7	6 3/8	162	17 3/8	441.3
13	6 ½	165.1	6 5/8	168.3	18	457.2
14	6 ¾	171.5	6 7/8	174.6	18 3/4	476.2

Sport Training



Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 ½	114.3	12 3/8	314.9
6	4 ¾	120.6	4 ¾	120.5	13 1/16	332.4
7	5	127	5	127	13 3/4	350
8	5 ¼	133.3	5 ¼	133.3	14 7/16	367.2
9	5 ½	139.7	5 ½	139.7	15 1/8	384.8
10	5 ¾	146	5 ¾	146	15 13/16	402.3
11	6	152.4	6 1/16	154	16 1/2	419
12	6 ¼	158.7	6 5/16	160	17 1/4	438.1
13	6 ½	165.1	6 9/16	166.7	17 7/8	454
14	6 ¾	171.5	6 13/16	173	18 5/8	473

Reverse Hind



Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 ¾	120.5	13	330
6	4 ¾	120.6	5 1/16	128.5	13 3/4	348.4
7	5	127	5 ¼	133.3	14 7/16	366.7
8	5 ¼	133.3	5 5/8	142.8	15 1/8	385
9	5 ½	139.7	5 7/8	149.2	15 7/8	403
10	5 ¾	146	6 1/8	155.5	16 9/16	421.6
11	6	152.4	6 7/16	163.5	17 3/8	441.3
12	6 ¼	158.7	6 11/16	170	18 1/8	460
13	6 ½	165.1	7	177.8	18 7/8	479.4
14	6 ¾	171.5	7 1/4	184.1	19 1/2	495.3

Suspensory Branch Front



Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 7/8	123.8	13 1/8	333.5
6	4 ¾	120.6	5 1/4	133.4	13 7/8	351.8
7	5	127	5 ½	139.7	14 9/16	370.3
8	5 ¼	133.3	5 3/4	146.5	15 5/16	389.1
9	5 ½	139.7	6 1/16	153.7	16	407.4

10	5 ³ / ₄	146	6 ⁵ / ₁₆	160.5	16 ³ / ₄	426
11	6	152.4	6 ⁵ / ₈	168.3	17 ¹ / ₂	444.5
12	6 ¹ / ₄	158.7	6 ⁷ / ₈	174.6	18 ¹ / ₄	463.5
13	6 ¹ / ₂	165.1	7 ¹ / ₈	181	19	482.6
14	6 ³ / ₄	171.5	7 ⁷ / ₁₆	190	19 ¹¹ / ₁₆	500

Bevel Lateral Roller




Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 7/8	123.8	13 1/4	336.6
6	4 ¾	120.6	5 1/8	130.2	14	355.6
7	5	127	5 7/16	138.1	14 11/16	373.1
8	5 ¼	133.3	5 11/16	144.5	15 7/16	392.1
9	5 ½	139.7	6	152.4	16 3/16	411.2
10	5 ¾	146	6 ¼	158.8	16 7/8	428.6
11	6	152.4	6 ½	165.1	17 5/8	447.7
12	6 ¼	158.7	6 13/16	173	18 3/8	466.7
13	6 ½	165.1	7 1/16	179.4	19 1/8	485.8
14	6 ¾	171.5	7 5/16	185.7	19 13/16	503.2

Suspensorion



Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	5	127	12 7/8	327
6	4 ¾	120.6	5 ¼	133.4	13 ½	342.9
7	5	127	5 ½	139.7	14 ¼	362
8	5 ¼	133.3	5 ¾	146.1	15	381
9	5 ½	139.7	6	152.4	15 ¾	400
10	5 ¾	146	6 5/16	160.3	16 7/16	417.5
11	6	152.4	6 5/8	168.3	17 1/8	435
12	6 ¼	158.7	6 7/8	174.6	17 7/8	454
13	6 ½	165.1	7 1/8	180.9	18 9/16	471.5
14	6 ¾	171.5	7 7/16	188.9	19 ¼	489

Performance



Size	Width	In Millimeters	Toe to Heel Height	In Millimeters	Circumference	In Millimeters
5	4 ½	114.3	4 ½	114.3	12 3/8	314.9
6	4 ¾	120.6	4 ¾	120.5	13 1/16	332.4
7	5	127	5	127	13 3/4	350

8	5 ¼	133.3	5 ¼	133.3	14 7/16	367.2
9	5 ½	139.7	5 ½	139.7	15 1/8	384.8
10	5 ¾	146	5 ¾	146	15 13/16	402.3
11	6	152.4	6 1/16	154	16 1/2	419
12	6 ¼	158.7	6 5/16	160	17 1/4	438.1
13	6 ½	165.1	6 9/16	166.7	17 7/8	454
14	6 ¾	171.5	6 13/16	173	18 5/8	473