

Molift Mover 205 Low Base

Molift Mover 205 is a flexible all-round lifter for all lifting and moving situations in the home, care institutions and hospitals. It lifts from the floor without adjusting the height of the column and has an excellent lifting height. Its lightweight, low chassis and wheel positioning make the Molift Mover 205 easy to maneuver.



Properties

- Easy to maneuver
- Excellent lifting height
- Fits easily under furniture and beds
- 4-point suspension with rapid connection system
- Environmentally friendly batteries
- Height of legs, only 2 3/8 inches

LIFTABILITY

molift[®]
by Etac

Molift Mover 205 Low Base

Molift Mover 205 is supplied with 4-point suspension as standard. Used together with a Molift Sling, it provides comfortable and correct positioning from lying to sitting and vice versa. Molift Mover 205 can also be supplied with optional suspension units which can be changed quickly and easily with a rapid connection system. Molift Mover 205 is equipped with automatic service notification.

This consists of integrated software which calculates when the lifter requires a service. It is also equipped with electronic and manual emergency lowering and emergency stop. Molift Mover 205 is available with an extra low chassis, only 6 cm.



Technical facts

Item no.

xxxx

Lifting capacity

205 kg

Total weight

39.3 kg

Lifting height-Min 340 mm –
Max 1800 mm**Chassis Width**

700 mm

Leg height

115 mm / 60 mm low chassis

Battery

NiMH 14.4 V - 2.2 Ah

Molift Mover 205 can be used with all passive Molift slings designed for 4-point suspension.

Accessories

**Molift Mover 205 Scale**

The scale complies with all applicable technical standard and is approved class III scale. The scale can calculate Body Mass Index.

Item no: 1431010

**Molift Stretcher**

A scoop carrier specially adapted to users with back injuries. It is X-ray compatible and can be suspended directly from the standard 4-point suspension.

Item no: 2150103

**Molift Mover Walking Arms**

Molift Mover is available with adjustable walking arms. The walking arms can be used together with the Molift ambulating vest for walking exercises. The arms fold along the column when not in use.

Item No.: 1411002