

molift QuickRaiser 205



Molift QuickRaiser 205 is a state of the art sit-to-stand hoist, with an **electrically adjustable leg base** and an **outstanding lifting capacity**. The standing hoist is comfortable, stable and safe with **exceptional manoeuvrability**. The unique wheel position, angled push bar and the extraordinary weight distribution of the base provides a **small turning radius** and easy manoeuverability even in very narrow spaces. With the small castor option (Ø 50mm), it is only 7,5 cm high (3"), and fits under most beds and chairs while the electrically adjustable legs enables accessibility almost everywhere.

The Active Lifting Arm (four point spreader bar) and the RgoSling Active provides **extra support** when hoisting from sit to stand. The inclined column **encourages a natural pattern** of movement and fulfills the hoisting needs of users with some weight bearing capacity.

LIFTABILITY









molifti QuickRaiser 205

Features

Ergonomic Design Slanting column lifts the user up and forwards to imitate the natural movement pattern.

Comfort, Safety and Stability Active Lifting Arm (four point spreader bar) and the RgoSling Active provides extra support when hoisting from sit to stand.

Automatic Warning Annual periodic inspection/service.

Easy to Manoeuver Easy to operate in moving and handling situations; toilet, sit-to-stand, stand-up training.

NEW Operation Panel and hand control provides intuitive and safe operation. Also includes a quick guide for sling application and safe use.

Technical Facts

Item No.

29000 (Active 50/75mm) 29001 (Standup 50/75mm) 29100 (Active 75/100mm) 29101 (Standup 75/100mm) Max User Weight 205kg Total Weight 44kg Chassis 17kg Lifting Arm 4.5kg Leg Support 3.5kg

Chassis 1/kg Lifting Arm 4.5kg Leg Support 3.5kg Column (heaviest part) 19kg

Battery SLA 24V **Charge Time** 4hours

No. of Lifts (75kg) with Fully Charged Battery >100 Lifting Speed 3.5cm/s 75kg 3.0cm/s 205kg

Lifting Range 60cm Width of Wheelbase - External 62cm - Internal 20cm Height of Wheelbase Castor Ø 75/50 76mm

Height of Wheelbase Castor Ø 100/75 101mm

Ground Clearance Castor Ø 75/50 27mm Ground Clearance Castor Ø 100/75 52mm

Height of Footplate Castor Ø 75/50 65mm

Height of Footplate Castor Ø 100/75 90mm

Turning Radius 115cm

Cleaning Instructions Use a solvent-free cleaning agent (ph Level 5-9) or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F

Accessories

Stand Up Lifting Arm 2910061
Ideal to use with RgoSling Standup slings.

Adjustable Leg Support 2910055
The leg support is adjustable without tools.

Inlays for Footplate 2920151 set of 2
Inlays with slip resistant material to increase friction.

