

# FOAM HEEL DRESSINGS

Vitále® Foam Heel Dressings are meticulously designed to provide superior protection and comfort for heel wounds. These dressings feature a highly absorbent foam construction that conforms easily to the unique shape of the heel, ensuring full coverage and protection. The foam not only absorbs exudate effectively but also retains warmth, promoting an optimal healing environment. The vapor-permeable film backing allows air to circulate while preventing moisture and external contaminants from entering, thereby reducing the risk of maceration and infection. The unique design of Vitále® Foam Heel Dressings allows them to be easily applied and securely held in place, with overlapping sides that help seal off any potential leaks. These dressings are ideal for treating wounds with moderate to heavy exudate, providing reliable care for patients in both clinical and home settings.



## USAGE

Type: Primary

Stages: II-IV

Drainage/Exudate: Moderate to Heavy

### Wound Types:

Pressure Sores

Venous Ulcers

Diabetic Ulcers

Grafts & Skin Donor Sites

Surgical/Post-Op

Traumatic Wounds

Burns

## FEATURES & BENEFITS

- **High Absorbency:** The foam construction efficiently absorbs and retains exudate, reducing the risk of maceration and promoting a moist healing environment.
- **Conformable Design:** The dressing is designed to conform to the unique shape of the heel, ensuring comprehensive coverage and protection.
- **Vapor-Permeable Film:** The breathable film backing allows moisture vapor to escape while protecting the wound from external contaminants.
- **Secure Fit:** The dressing's overlapping sides ensure a secure fit, sealing off potential leaks and providing consistent protection throughout the wear time.
- **Versatile Application:** Suitable for a variety of heel and ankle wounds, including pressure sores, arterial ulcers, diabetic ulcers, post-operative wounds, and second-degree burns.

Vitále® Foam Heel Dressings are trusted by healthcare professionals for their advanced design and reliable performance. The high absorbency, conformable fit, and breathable protection ensure that heel wounds are managed effectively while maintaining patient comfort. These dressings are versatile, easy to apply, and ideal for a wide range of heel and ankle wounds, making them an essential tool in both clinical and home care settings. By choosing Vitále® Foam Heel Dressings, healthcare providers can offer superior wound care that enhances healing and improves patient outcomes.

## ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION	SIZE	UNITS/BOX	HCPC CODE
20400	Vitále® Foam Heel Dressing	7" x 7" (18 cm x 18 cm)	5	A6210