


VITÁLE®

NON-ADHESIVE FOAM DRESSINGS

Vitále® Non-Adhesive Foam Dressings are expertly crafted to provide optimal wound care by maintaining a moist healing environment and offering high absorbency. These dressings feature a two-layer construction: an absorbent foam layer that retains exudate and a waterproof yet vapor-permeable film layer (except for the 20463 model) that protects the wound while allowing moisture vapor to escape. This design ensures effective exudate management, preventing bacterial contamination and skin maceration. The soft and flexible nature of these dressings makes them comfortable on fragile skin, and they can be worn for up to seven days. Vitále® Non-Adhesive Foam Dressings are ideal for a variety of wound types with moderate to heavy exudate, making them a reliable choice for both clinical and home care settings.



USAGE

Type: Primary

Stages: II-IV

Drainage/Exudate: Moderate to Heavy

Wound Types:

Pressure Sores

Venous Ulcers

Diabetic Ulcers

Grafts & Skin Donor Sites

Surgical/Post-Op

Traumatic Wounds

FEATURES & BENEFITS

- **Highly Absorbent:** The foam layer efficiently absorbs and retains exudate, maintaining an optimal moist wound environment.
- **Moisture and Bacteria Barrier:** The waterproof, vapor-permeable film (on applicable models) protects the wound from external contaminants while allowing excess moisture vapor to escape, reducing the risk of maceration.
- **Soft and Flexible:** The dressings are designed to be gentle on fragile skin, providing a comfortable fit and reducing trauma during dressing changes.
- **Long Wear Time:** Can be worn for up to seven days, reducing the frequency of dressing changes and enhancing convenience for both patients and caregivers.
- **Non-Adhesive:** Ideal for patients with sensitive skin or those who are allergic to adhesives, ensuring a pain-free removal process.

Healthcare professionals trust Vitále® Non-Adhesive Foam Dressings for their exceptional ability to manage exudate and protect wounds. The advanced two-layer design ensures that wounds are kept in an optimal moist environment, promoting faster healing and reducing the risk of infection. Clinical feedback highlights significant improvements in healing times and patient comfort, with the non-adhesive design minimizing irritation and providing secure protection. These dressings are versatile, easy to apply, and ideal for a wide range of wound types, making them a preferred choice in both clinical and home care settings. Vitále® Non-Adhesive Foam Dressings are an essential component of advanced wound care, leading to better patient outcomes and satisfaction.

ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION	SIZE	UNITS/BOX	HCPC CODE
20463	Vitále® Non-Adhesive Foam Dressing	3"x 3" (7.5 cm x 7.5 cm)	10	A6212
20444	Vitále® Non-Adhesive Foam Dressing	4" x 4" (10 cm x 10 cm)	10	A6212

