

GOT MORE QUESTIONS?

Check out our list of Frequently Asked Questions
at vhealth.link/c647e for helpful answers.



And if that doesn't answer your question, our customer
service team would love to help! Feel free to connect with
them by phone, e-mail, or chat on our website.



service@vivehealth.com



1-800-487-3808



vivehealth.com

vive®



TRANSFER BELT WITH HANDLES

Owner's Manual
RHB1011N

vivehealth.com

OVERVIEW


We are constantly answering questions and recording helpful videos to make using your Vive Transfer Belt as easy as possible. Check out the included link and QR code to help you through the process.




To see all of the FAQ's in one place visit vhealth.link/c647e

QUICK START GUIDE

1. Insert strap into the buckle through the end with the teeth.



2. Slide strap through the hole opposite the teeth.



QUICK RELEASE

To quickly release the strap, pull the release tab and slide the strap free.

