

## DropBot Chip Product Insert

Thank you for purchasing one or more verified Sci-Bots DropBot Chips. Follow the instructions given in your DropBot Instruction Manual to load and operate the chip.

Each chip has been carefully tested to ensure that it passes our rigorous quality control procedures. A nominal “driving voltage” is indicated on the packaging – we guarantee that at this voltage, droplets of our test solution (PBS + 0.1% T90R4) should move over all of the electrodes on the array.

Contact us at [support@sci-bots.com](mailto:support@sci-bots.com) with any questions.

