

ELMI CM-7S Plus

Swing-Out Centrifuge



FEATURES



Cold room and refrigerator safe



Compact, Lightweight, and Versatile



Rotor imbalance technology



Automatic lid lock for convenience and safety



Rapid acceleration/braking system



Multiple brake speeds, 30 seconds to 2.5 minutes



Digital timer



LCD Control Panel

Control speed in RPM or G-Force

SPECIFICATIONS

| | |
|--|------------------------------------|
| Rotor Rotation Speed | 100 - 3,500 RPM |
| Timer | 1 second - 12 hours |
| Maximum Centrifugal Force | 2,300 RCF (G-Force) |
| Brake Levels | 6 |
| Working Temperature Range | from +10 to +40 °C |
| Relative Air Moisture at 20 °C (68 °F) | 80% |
| Voltage and Frequency | 100 - 240 V; 50/60 Hz |
| Maximal Consumable Power | 250 W |
| Size (L x W x H) | 431 x 410 x 247 mm 17 x 16x 10" |
| Weight | 13.2 kg/29 lbs |

ACCESSORIES

| SKU | Product Name | PRICE |
|-----------|---|---------|
| CM-6MT-8 | CM-7S PLUS Benchtop Swing-Out Centrifuge, 3500 RPM. Stainless Steel | \$1,450 |
| CM-6MT-3S | Rotor 6M.05 - 12 x 15ml Test Tubes - Steel adapters included | \$254 |
| CM-6MT-5S | Rotor 6M.02 - 24 x 15ml Test Tubes - Steel adapters included | \$362 |
| CM-6MT-7 | Spacer for microtainers (3.5-10ml)-Compatible with 6M.05 & 6M.02 | \$3 |
| CM-6MT-20 | Rubber spacer for glass tubes-Compatible with 6M.05 & 6M.02 | \$2 |
| CM-6MT-4 | Rotor 6M.01 - 4 x 50ml Test Tubes - Plastic adapters included | \$190 |
| CM-6MT-12 | Rotor 6M.06 - 6 x 50ml Test Tubes - Plastic adapters included | \$289 |
| CM-6MT-15 | Adapter/Insert for 15ml tubes - for use with 50ml buckets from rotors 6M.01&6M.06 | \$15 |
| CM-6MT-21 | Combo Rotor 6M.07 2x50ml + 6x15ml - Adapters Included | \$299 |
| CM-6MT-22 | Rotor 6M.08 4 x 60ml Test Tubes - Stainless steel adapters included | \$295 |

Distributed By:



medicalspsasupply.com
(800) 880-8260