






























Glossary of symbols

Self-adhesive electrodes

	Manufacturer		Non-sterile
	Follow instructions for use		Medical device
	Minimum shelf life		Batch code
	Attention		Article number
	Free for sale		Unique product identification number
18+	Only for people over 18 years of age		Temperature limit
	Protect from moisture		Protect from sunlight
	CE marking		Use only on healthy skin
	Do not use by people with pacemakers or defibrillators		Do not use on the head, carotid artery or heart
	Do not apply over implants		Do not use while driving a vehicle
	May only be used by one person		Never stick electrodes on top of each other
	Keep away from children		Amperage (I) not higher than 0.25 W / cm²
	Humidity restriction		For plug-in connections with a receptacle of Ø 1.9 - 2.1 mm
	Free of latex		For push-button connections with 2.2mm, 3.2mm, 3.5mm or 3.8mm
	Read the instructions for use of the electrical stimulation device before use		Electrodes for use on the skin with a TENS / EMS device