



DIAGNOSTIC & OPERATIVE SHEATH



Hysteroscopic Examination Sheath With Atraumatic Distal End & Rotatable Stop Cock For 4mm & 2.9mm 30' Scope



Hysteroscopic Operative Sheath, Contributions With 2 Detachable Suction Irrigation Plugs For 4mm & 2.9mm 30' Scope



WORKING ELEMENT MONOPOLAR & BIPOLAR

- · Monopolar Resectoscope Can Be Used In Complex Intrauterine Operation Without Any Difficulty.
- · Bipolar Resectoscope Are Particularly Well Sulted For Infertility Patients.
- · Bipolar Resection Benefits From Generation Of Plasma To Achieve Outstanding Cutting Power.
- By Using Hf Current Several Operations Like Hysteroscopic Endometrium Ablation,
- Myoma Resection, Adhesiolysis Or Septum Dissection Can Be Performed.
- · An Electrolyte Free Irrigation Fluid Is Used As Distension Media

WORKING ELEMENT





Outer Sheath



Inner sheath



Working Element Passive Monopolar / Bipolar



Obturator

ELECTRODES



Electrode Roller-Ball



Electrode Cutting Loop,



Electrodes Knife

www.incisioninstrument.com

incisioninstrument@gmail.com

+91 9892736921

shop no 2 , malik compound , opp white city by lodha , lokhandwala , akurli road , kandivali (e) , mumbai 400101 s.