

## Type IIR disposable medical mask

High splash resistant mask for enhanced protection



Manufactured to European standards (EN14683:2019 +AC:2019) these type IIR masks are designed for protecting people from infectious diseases.

Conforming to the highest standards of protection, they are suitable for use by the NHS and recommended by World Health Organisation for the healthcare sector. The 'R' specifies that they are fluid resistant to offer increased protection compared to a standard mask as they block potentially-infected droplets from entering the respiratory system.

## Available as one size fits all and supplied in a box of 50 pieces for single use.

- **©** CE marked
- ☑ Lightweight and breathable 3-ply protection
- High splash resistance
- **⊘** Secure ear loop fit
- Omfortable integrated nose clip
- **⊙** High Bacterial Filtration Efficiency (BFE) ≥ 98%

## How to use and discard of a disposable mask:



Clean hands with soap and water or an alcohol hand sanitiser



Cover mouth and nose with the mask – make sure there are no gaps



Avoid touching the mask, if you do wash your hands or hand sanitise



Replace the mask as soon as it gets damp - do not re-use single use masks



Remove the mask from behind – do not touch the front



Discard into a closed bin



Clean hands with soap and water or hand sanitiser