



More information on the website
radwag.com/en/info,w1,KRZ

SAL/C Laboratory Antivibration Table



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	72.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MUYA Modules
MAS1 Modules
MA 5Y Moisture Analyzer

MPS Modules
C32.PM Precision Balance
MA X7.A Moisture Analyzer
MA 5Y.IC Moisture Analyzer

Device dimensions

SAL/C

