



Instructions for Use

AA00003-IFU

Rev. Date: Sep-16

Revision: 1

Page 1 of 1

Monoclonal Mouse Antibody to Human B Cell (CD20)

Description:

Immunogen:	BALB/C mice were injected with human tonsil B cells.
Clone:	L26.
Isotype:	IgG2a, Kappa.
Format:	This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further dilution is required.
Specificity:	This antibody is specific to a 33 kD polypeptide present on the majority of B cells in peripheral blood and lymphoid tissue and also with a minor component of 30 kD. No reactivity with other haematopoietic cells has been observed.

Uses/Limitations:

Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage:

2-8° Centigrade.

Procedure:

We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions:

Contains Sodium azide as a preservative.
Do not pipette by mouth.