

A50 Cooper II

US Patent 505062
Suggested LCode: L1906

APPLICATIONS:

The Cooper II Ankle Stabilizer brace offers bilateral support and can be worn on either foot, which eliminates the need for separate left and right versions and simplifies inventory management.

To better control inversion and eversion and other instabilities of the ankle, the medial and lateral cinch straps originate from the forefoot. Additionally, this stabilizer is both comfortable to wear and easy to apply, making it an ideal solution for any ankle or forefoot needs. From emergency room treatment to the field of competition, the Cooper II Ankle Stabilizer offers reliable support for a seamless transition.

FEATURES:

- Very Comfortable
- Breathable, perforated material
- Bilateral support (can be worn on left or right foot)
- Non-stretch medial and lateral straps help prevent rollover
- Easy donning and doffing

INDICATIONS:

- Inversion or eversion ankle sprains
- Medial or lateral ankle instability
- Long-term management to help reduce chance of frequent ankle injuries
- Subtalar joint instability
- Syndesmotic ankle sprain (high ankle sprain)

Washing Instructions:

Hand wash in cold water with mild detergent.
NO chlorine bleach. Air dry. Do NOT put in dryer.



ANKLE SUPPORTS (By men’s shoe size)
Size may vary according to ankle circumference.

X-Small	5 - 6
Small	6 - 7
Medium	7 - 9
Large	10 - 12
X-Large	12 - 14
2X-Large	14 - 16

1850 Diplomat Drive, Suite #100
Farmers Branch, TX 75234
PH#: 214-638-6422 FX#: 214-638-6425

