

E10-U "Gel*Bow"

Universal Elbow Strap with Gel Pad
Suggested LCode: L3999

FEATURES:

- Circumferential loop and lock counterforce strap with tongue
- Covered pressure gel pad to alleviate stress over the medial/lateral epicondyle
- Gel pad can be chilled in a freezer at night or during the day and applied as cold therapy as needed.
- 1/8" latex free open cell foam with one side UBL, one side nylon
- 1-1/2" strap width, 3" bolster gel pad
- Universal sizing, adjustable to fit 7"-18" forearm circumference

INDICATIONS:

- Sprains/Strains
- Medial/Lateral Epicondylitis
- Generalized Tendomyopathy

CONTRAINDICATIONS:

DO NOT wear if you have circulatory problems or are susceptible to dermatitis. If rash, blisters, or skin irritations occur, consult your physician.



Washing Instructions:

Hand wash in cold water with mild detergent.
NO chlorine bleach. Air dry. Do NOT put in dryer.



1850 Diplomat Drive, Suite #100
Farmers Branch, TX 75234
PH#: 214-638-6422 FX#: 214-638-6425

