



# INSTRUCTIONS FOR USE

## CALF SLEEVE, C1



**1.** Position support with wide end up and logo forward, seam in back. With thumbs inside of the wide end, grasping outside with fingers, pull the support over foot and up calf until it feels snug. For ease of application, try to pull with one continuous motion.



**2.** For best fit, position the center seam along the center of the calf.  
FITTING TIP: depending on the length of the calf versus the length of the support, the top of the support can be folded down and “cuffed” for a more snug fit.

Still need help?  
Scan here to watch a video tutorial!



C1