

K42-HT "Rehab Knee"

Suggested LCode: L1832

FEATURES:

- Anterior closure for easy donning and doffing
- Patellar opening for positioning and comfort
- · Lightweight, 1/8" nylon 2-sides neoprene
- · Open popliteal area
 - Avoids binding and pinching for sensitive patients
 - Promotes supplemental air contact for cooling and comfort
- Two upper and lower circumferential straps
- DJ ROM hinge for locking out or set for flexion and extension control
- · 13" in length
- · Hinge arms made of aluminum

INDICATIONS:

- Post-operative management and rehabilitation
- Medial and lateral collateral ligament instability
- · Preventative care



Washing Instructions:

Hand wash in cold water with mild detergent. NO chlorine bleach. Air dry. Do NOT put in dryer.



SIZING

Measure knee circumference with leg extended.

>	K-Small	12" - 13"
5	Small	13" - 14"
١	Medium	14" - 15"
L	_arge	15" - 16"
>	X-Large	16" - 18"
2	2X-Large	18" - 20"
3	3X-Large	20" - 22"
4	4X-Large	22" - 24"
(Custom	24"+ UP

1850 Diplomat Drive, Suite #100 Farmers Branch, TX 75234

PH#: 214-638-6422 FX#: 214-638-6425

