

LC10 Lumbosacral Corset Orthosis

PDAC Approved A4466

FEATURES:

- KoolFlex™ posterior panel
- Latex Free, Elastic anterior panels with hook and loop closure
- Four removable pre-contoured rigid steel posterior stays
- Two reinforced circumferential straps
- Heavy duty reinforced loops
- Anterior hook and loop Velcro® closure
- Posterior height of 10"

INDICATIONS:

- Lower back pain
- Post-operative care
- Muscle Spasms
- Facet Syndrome

FITTING NOTES:

- When fitting the corset, first secure the hook & loop closure on the body.
- Be sure that all straps are pulled straight and evenly through the loops.
- Fasten the hook & loop straps and adjust the buckles to the desired level of compression.



Washing Instructions:

Hand wash in cold water with mild detergent.
NO chlorine bleach. Air dry. Do NOT put in dryer.

Front



Back



SIZING

Measure around waist.

X-Small	25" - 32"
Small	28" - 35"
Medium	32" - 41"
Large	34" - 44"
X-Large	40" - 48"
2x-Large	42" - 54"
3x-Large	44" - 59"
Custom	50" & UP

1850 Diplomat Drive, Suite #100
Farmers Branch, TX 75234
PH#: 214-638-6422 FX#: 214-638-6425

