

KnitRite AFO Interface Socks

To select the correct size, measure your leg from the tip of the toe to the top of the AFO / KAFO. Select the sock which is closest to that length.

There is no heel pocket, so the shoe size doesn't matter. The sock is like a long stretchy tube that will mould to the contours of the foot and leg. Measurements are the unstretched socks.

SmartKnit AFO Interface Seamless Socks



Size	Width	Length
Infant Small AFO	3.5cm	18cm
Infant AFO	5cm	23cm
Child Small AFO	5.5cm	28cm
Child AFO	6.5cm	36cm
Child Extra Long AFO	7.5cm	51cm
Adult Small AFO	7.5cm	46cm
Adult Regular AFO	9cm	51cm

Patterned AFO Core-Spun Interface Socks



Size	Width	Length
Child Small	5.5cm	28cm
Child Regular	6.5cm	36cm
Child Extra-Long	7cm	51cm
Adult Small	7.5cm	46cm
Adult Regular	9cm	51cm



KnitRite AFO & KAFO Interface Socks

To select the correct size, measure your leg from the tip of the toe to the top of the AFO / KAFO. Select the sock which is closest to that length.

There is no heel pocket, so the shoe size doesn't matter. The sock is like a long stretchy tube that will mould to the contours of the foot and leg. Measurements are the unstretched socks.

SmartKnit Big Toe AFO Seamless Socks



Size	Width	Length
Infant	5cm	23cm
Child Small	5.5cm	28cm
Child Regular	6.5cm	36cm
Adult Small	8cm	36cm
Adult Regular	9cm	51cm

SmartKnit KAFO Interface Seamless Socks



Size	Width	Length
Child XX-Long KAFO	6.5cm	46cm
Adult Narrow KAFO	10 x 6.5cm	61cm
Adult Regular KAFO	15 x 9cm	71cm



Torso Interfaces

Outstanding stretch of over 150% allows the fit to be virtually wrinkle-free.

Style options include: **axilla flaps** that protect delicate underarm skin from rubbing and abrasion. Another style option is the weight of the garment. The **original thickness** provides additional protective comfort, while **lightweight thickness** is cooler to wear.

Torso Interface® Turtleneck with Sleeves Size Width (unstretched) **Fits Circumference** Length X-Small 10cm 30cm 25-30cm 35cm 35-45cm Small 15cm Medium 20cm 40cm 45-60cm

Torso Interface® Without Straps					
	Size	Width (unstretched)	Length	Fits Circumference	
	Small	15cm	20cm	35-45cm	
	Medium	20cm	40cm	45-60cm	
	Medium Long	20cm	50cm	45-60cm	
	Large	25cm	50cm	55-75cm	
	Large Long	25cm	60cm	55-75cm	
	X-Large	30cm	60cm	66-95cm	
	KS	35cm	60cm	76-122cm	

Torso Interface® Crew Neck: With Sleeves					
	Size	Width (unstretched)	Length	Fits Circumference	
	Small	15cm	35cm	35-45cm	
	Medium	20cm	40cm	45-60cm	
	Medium Long	20cm	45cm	45-60cm	
	Large	25cm	50cm	55-75cm	
	Large Long	25cm	60cm	55-75cm	
	X-Large	30cm	60cm	66-95cm	
	X-Large Long	30cm	71cm	66-95cm	



Standard Torso Interfaces

Outstanding stretch of over 150% allows the fit to be virtually wrinkle-free.

Style options include: **axilla flaps** that protect delicate underarm skin from rubbing and abrasion. Another style option is the weight of the garment. The **original thickness** provides additional protective comfort, while **lightweight thickness** is cooler to wear.

Torso Interface® Crew Neck



Size	Width (unstretched)	Length	Fits Circumference
X-Small	10cm	30cm	25-30cm
Small	15cm	35cm	35-45cm
Medium	20cm	40cm	45-60cm
Medium Long	20cm	56cm	45-60cm
Large	25cm	50cm	55-75cm
Large Long	25cm	60cm	55-75cm
X-Large	30cm	60cm	66-95cm
KS	35cm	60cm	76-122cm
KX	35cm	71cm	76-122cm

Torso Interface® V Neck



Size	Width (unstretched)	Length	Fits Circumference
Small	15cm	35cm	35-45cm
Medium	20cm	40cm	45-60cm
Medium Long	20cm	45cm	45-60cm
Large	25cm	50cm	55-75cm
Large Long	25cm	60cm	55-75cm
X-Large	30cm	60cm	66-95cm
KS	35cm	60cm	76-122cm

Torso Interface® Lightweight V Neck



Size	Width (unstretched)	Length	Fits Circumference
Small	15cm	40cm	41-69cm
Medium	20cm	45cm	56-94cm
Medium Long	20cm	56cm	56-94cm
Large	25cm	50cm	76-127cm
Large Long	25cm	60cm	76-127cm
X-Large	30cm	60cm	92-153cm

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.store	09477873	help@sensorysmart.store	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands, B73 5TE, UK



Left Axilla Torso Interfaces

Outstanding stretch of over 150% allows the fit to be virtually wrinkle-free.

Style options include: **axilla flaps** that protect delicate underarm skin from rubbing and abrasion. Another style option is the weight of the garment. The **original thickness** provides additional protective comfort, while **lightweight thickness** is cooler to wear.

Torso Interface ® V Neck: Left Axilla Flap Size Width (unstretched) Length **Fits Circumference** Small 15cm 35cm 35-45cm Medium 20cm 40cm 45-60cm 25cm 50cm 55-75cm Large Large Long 25cm 60cm 55-75cm 60cm 66-95cm X-Large 30cm

Torso Interface® Lightweight V Neck: Left Axilla Flap					
	Size	Width (unstretched)	Length	Fits Circumference	
	Small	15cm	40cm	41-69cm	
	Medium	20cm	45cm	56-94cm	
	Medium Long	20cm	56cm	56-94cm	
4 3	Large	25cm	50cm	76-127cm	
15 E	Large Long	25cm	60cm	76-127cm	
	X-Large	30cm	60cm	92-153cm	



Double Axilla Flap Torso Interfaces

Outstanding stretch of over 150% allows the fit to be virtually wrinkle-free.

Style options include: **axilla flaps** that protect delicate underarm skin from rubbing and abrasion. Another style option is the weight of the garment. The **original thickness** provides additional protective comfort, while **lightweight thickness** is cooler to wear.

Torso Interface® V Neck: Double Axilla Flaps Size Width (unstretched) Length **Fits Circumference** Small 35-45cm 15cm 35cm Medium 20cm 40cm 45-60cm Medium Long 20cm 45cm 45-60cm 25cm 50cm 55-75cm Large 55-75cm Large Long 25cm 60cm X-Large 30cm 60cm 76-122cm

Torso Interface® Lightweight V Neck: Double Axilla Flaps					
	Size	Width (unstretched)	Length	Fits Circumference	
	Small	15cm	40cm	41-69cm	
	Medium	20cm	45cm	56-94cm	
	Medium Long	20cm	56cm	56-94cm	
	Large	25cm	50cm	76-127cm	
E =	Large Long	25cm	60cm	76-127cm	
	X-Large	30cm	60cm	92-153cm	

Torso Interface® Crew Neck: Double Axilla Flaps						
	Size	Width (unstretched)	Length	Fits Circumference		
	Small	15cm	35cm	35-45cm		
	Medium	20cm	40cm	45-60cm		
	Medium Long	20cm	56cm	45-60cm		
	Large	25cm	50cm	55-75cm		
	Large Long	25cm	60cm	55-75cm		
	X-Large	30cm	60cm	66-95cm		

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.store	09477873	help@sensorysmart.store	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands, B73 5TF, UK