

Spio Size Chart


SPIO TLSO System (Vest)

Choose size based on chest measurement.
 Take measurement with patient lying down.

Available Sizes:

V 41 V 44 V 47 V 51 V 54
 V 60 V 67 V 70 V 73 V 77

Chest
 circumference in
 centimetres at
 nipple line




SPIO Upper Body Orthosis (Shirt)

Choose size based on chest measurement.
 Take measurement with patient lying down.
 Short and long sleeve available.

Available Sizes:

S 41 S 44 S 47 S 51 S 54 S 57
 S 60 S 64 S 67 S 70 S 73 S 77

Chest
 circumference in
 centimetres at
 nipple line




SPIO Lower Body Orthosis (Pant)

Choose size based on waist measurement.
 Take measurement with patient standing up.
 Short and long pant available.

Available Sizes:

P 41 P 43 P 45 P 47 P 49 P 51
 P 53 P 55 P 57 P 59 P 61 P 63
 P 65 P 67 P 69 P 71 P 73

Waist
 circumference in
 centimetres at
 navel



SPIO Wrist Hand Orthosis (Glove)

Choose size based on palm measurement.

Available Sizes:

G 4 G 4.5 G 5 G 5.5 G 6 G 6.5
 G 7 G 7.5 G 8

Measure across
 palm in centimetres

