


Product Brochure

Made for the modern dentist



 firststar






Reliable. Purposefully designed. Incredible value.

We got tired of choosing only two so we decided to fix that. For over 20 years, we've been designing and building dental equipment for every dentist, whether you're after the top of the line with a feature list that gives a CVS receipt a run for its money, or something a bit more budget friendly that still checks all of the boxes.

We set out to make dental equipment easy because we know you want to focus on what you do best - dentistry. With Firststar, it's not wrong to love your dental equipment, we won't tell anyone.





Features

- 1 Award-winning hydraulic system with soft start/stop.
- 2 New locking swing down armrests with single hand activation.
- 3 Dual articulating push button head rest for better access to the oral cavity.

FDC50 - The Evolution

The pinnacle in patient access, versatility, and comfort. The FDC50 is the big brain move when you're after the latest in patient comfort and practical ergonomics. We've been in the equipment game for some time, and the FDC50 is our pièce de résistance.

We pulled out all the stops when designing the 50, building upon our award winning hydraulic system adding soft start/stop and giving it the most modern recline movement.

Your patients will be lulled asleep and won't be awakened as you reposition them. Our new locking swing down armrests are disengaged with a single hand improving workflow and ergonomics.


Speaking of which, we also kept the dual articulating push button head rest further improving patient positioning to give you optimal access to the patient while you work. The FDC50 is so comfortable we wouldn't fault you for putting one of these in your living room, either as furniture or a piece of art.



Thin, tapered backrest gives ample room for your legs while working



Durable internal stops ensure consistent, reliable performance



Scan to find out more



FDC38 - Never go without

The FDC38 is a prime example of Firststar's build quality. You won't find plastics in any of the chair structure meaning maximum rigidity and longevity for you and your patients. Our classic thin aluminum backrest that tapers in width maximizes space for your legs while working the oral cavity but offers plenty of support for your patient.

With every FDC38, you and your practice get peace of mind, proven reliability, and incredible value. The only thing we don't include are the patients.





A solid workhorse that moonlights as a show pony, the FDC38 is feature rich and built upon Firststar's proven hydraulic drive system. Designed to maximize flexibility and patient access the FDC38 swivels 50 degrees and can be paired with either of our delivery systems or used on its own with a cabinet or cart mounted delivery.

Available in a wide selection of colors in both Naugasoft and optional Ultraleather®.



Incredible range of adjustability for the perfect fit



Scan to find out more



Pair our chairs with our Pivot Mount Delivery System


Our Pivot Mount Deliveries maximize workspace flexibility by allowing on the fly ambidextrous access to the patient and offer an extensive list of standard features.



Looking for a power optic system, integrated chair controls, and automatic handpiece control? Check. If you look, you'll also find thoughtful details such as needle bearings in the arms to eliminate play and ensure smooth operation for the lifetime of your product. The Pivot Mount Delivery System empowers you to do more with less.

View our available options to see how you can make our Pivot Mount Delivery even more powerful.

Scan to find out more



Our Pivot Mount Delivery System empowers better and better dentistry

... or one of our Post Mount Deliveries if that's more your flavor

If your practice calls for cuspidors and rinse cups we've got you covered. Our Post Mount Deliveries have all of the features in our Pivot Mounts with the addition of cuspidor controls bundled into a convenient package. With built in timed cup fill, bowl rinse, and the air and water hookups to supply your favorite accessory, our Post Mount Deliveries provide all the tools you need in an affordable package.

Supply your water how you want, either straight from the tap or from the standalone bottle system with the routing valve that's standard on all of our deliveries. Our Post Mount Delivery Systems put the power in your hands.





Built in timed cup fill, bowl rinse, and air and water QD's

Remote cuspidor controls (when equipped) integrated into Touch Pad

Scan to find out more



Take a closer look at our delivery systems

The industry leading value doesn't stop at the things you can see. We loaded our delivery systems with many "optional" features as standard, then built them out of premium parts and materials because at Firststar, it's the little things that matter.

Scan to find out more



All of our control heads feature:

- 1 Swing arm pneumatic brake
- 2 Handpiece coolant control
- 3 Automatic handpiece control
- 4 Integrated power optic system
- 5 Up to two 6 pin ISO-C and one four hole handpieces
- 6 Autoclavable three mode syringe (air, water, mist)
- 7 Tap/water bottle routing valve
- 8 Master power safety switch

Pair our deliveries with an assistant's instrumentation suite

- 1 Adds another three mode syringe
- 2 Includes 1x saliva ejector and 2x HVE's (1x HVE when equipped with cuspidor)
- 3 Integrated touch pad control
- 4 Replaceable solids collector
- 5 Compatible with optional hands free extra oral suction holder



Lights

Enhance your ability to assess and diagnose with our LED lights that help you see more clearly with truer colors. Our lights benefit your whole team by offering enhanced illumination of the oral cavity and energy efficiency over traditional halogen bulbs paired with ergonomic mounting and smooth three axis head movement.


Our RL-206 and RL-208 lights also feature a composite mode, giving you additional time when working with composites.



4 bulbs



6 bulbs



Improve your infection control with our LED lights, featuring touch free sensor activation and intensity adjustment. With the wave of a hand you can change the brightness level, switch between standard illumination and composite mode, and turn your light on and off. Providing an estimated 40,000 hours of energy efficient bulb life, our LED dental lights will illuminate your path to successful treatment.



8 bulbs

Scan to find out more



Doctor and Assistant Stools

Don't get so caught up thinking about your patients only to forget about your team! We offer plush, ergonomic doctor and assistant stools in the same colors as our dental chairs so you can complete the look.

Available in all of the same colors as our chairs in both Naugasoft (standard) and optional Ultraleather®, our stools will complement the look of your operator, and your team will thank you for it.



Scan to find out more



Chair	FDC50	FDC38
Quiet hydraulic lift and recline power system	X	X
Soft start/stop drive technology	X	
Chair swivel, 25° each side	X	X
Dynamic Pivot advanced recline movement cradles the patient and eliminates bunching at the lower back while reclining	X	
Full metal frame and base assemblies	X	X
Dual articulating, height adjustable head rest	Button activated	Button activated
Locking swing down arm rests	Single handed button lock mechanism	Lever lock mechanism
Standard upholstery	Ultraleather®	Naugasoft
Chair at lowest point	18"	18"
Maximum lift range	13"	13"
Max load rating	300 lbs (with delivery); 450 lbs without delivery	300 lbs (with delivery); 450 lbs without delivery


Delivery System	Pivot Mount	Post Mount
Compatible with FDC38	X	X
Compatible with FDC50	X	X
Four position handpiece mount with automatic control	Optional	Optional
Three position handpiece mount with automatic control	X	X
Ambidextrous access	X	
Three mode syringe (air, water, mist)	Standard	Standard
Power optic system (power supply, connectors, tubing)	Standard	Standard
12-24VAC power transformer for LED light and power optics	Included	Included
Self contained clean water system	X	X
Integrated chair controls	X	X
Swiveling cuspidor		Optional
Remote cuspidor controls		Included with cuspidor
Water routing valve	Standard	Standard
Air and water quick disconnects	Standard	Standard
Quick connect threaded handpiece tubing	2x Four hole and 1x 6 pin ISO-C (Up to 2x 6 pin ISO-C)	2x Four hole and 1x 6 pin ISO-C (Up to 2x 6 pin ISO-C)

Naugasoft Color Options

Choose from our most popular Naugasoft color options. For the complete list of available colors visit our website. Can't find the color you're looking for? Contact us today to see how we can help.



BLACK SATIN
PR 40




MULBERRY
PR 41



FAWN
PR 35



COCOA
PR 51



SADDLE
PR 52



DEEP SAPPHIRE
PR 55



BLUE FOG
PR 46



MINT LEAF
PR 45




Ultraleather® Color Options


Peak comfort starts with Ultraleather®. Choose from our popular colors below or view the full range on our website. Looking for something unique? Contact us today to learn more about our custom options.




PLUM
291-9385



BAYOU
291-2685




CYAN
291-2553




SONORA
291-4454



RAVEN WING
291-5607



GRANITE
291-5708



CHAMPAGNE
291-3719



STONE
291-3602



You focus on the dentistry. We'll handle the equipment.

Firststar may be a new name but we're no stranger to dental equipment. We've been manufacturing for some of the biggest names in the industry for over 20 years and now we're ready to fly solo. We know what makes the industry tick and also where it can improve which is why we're here to shake things up. Our goal is to simplify the equipment process - from selecting the right set up to getting parts and service. Regardless of which product you choose you can rest assured that we've got your back with thoughtful designs, proven reliability, and our standard warranty. And if you're ever in the neighborhood feel free to stop by and say hello, we don't bite!

18915 E Valley Highway
Kent, Washington 98032
United States

hello@firstardental.com
www.firstardental.com

