



FDC388

SIDE MOUNT PACKAGE

DELIVERY STANDARD FEATURES

- Traditional delivery system with automatic control for three handpieces
- Pre-installed power optic system
- Push button control arm brake
- Selectable water source (tap or water bottle)
- Swiveling cuspidor with timed bowl rinse and automatic cup filler
- Integrated heated water system with quick disconnect outlet
- Four bulb motion activated LED light with adjustable intensity and three axis light head movement
- Telescoping rear assistants suite with integrated controls



ADDITIONAL OPTIONS

- 4+1 handpiece configuration
- 3+1 rear assistant suite upgrade (2 HVE, 1 SE, 1 syringe)
- Hands free suction holder with extra HVE and tubing
- 6 or 8 bulb LED light and arm upgrade

Every hero needs a sidekick

FDC38 SIDE MOUNT PACKAGE

We combined our workhorse 38 chair with the extremely capable Side Mount Delivery System for those who need cuspidor functions. Like all of our deliveries, you get a plethora of standard features and our outstanding build quality with added water and cuspidor features. Like Batman's utility belt, the FDC38 Side Mount Delivery gives you all the right tools at your disposal.


CHAIR FEATURES

- Fully hydraulic drive system
- 50 degree chair swivel
- Synchronized seat and backrest movement
- Dual locking swing down armrests
- Thin, tapered aluminum backrest
- Four position programmable chair memory
- Dual articulating, height adjustable headrest
- Durable and stylish Naugasoft upholstery



FIRSTAR DENTAL CO.

18915 E Valley Hwy
Kent, WA 98032
hello@firstardental.com
www.firstardental.com



Reliable. Purposefully designed. Incredible value. We got tired of choosing only two so we decided to fix that. For over 20 years, we've been designing and building dental equipment for every dentist, whether you're after the top of the line or something a bit more budget friendly that still checks all the boxes. We set out to make dental equipment easy because we know you want to focus on what you do best - dentistry. With Firstar, it's not wrong to love your dental equipment.

