



# FDC50

PIVOT MOUNT PACKAGE

## DELIVERY STANDARD FEATURES

- Traditional delivery system with automatic control for three handpieces
- Pre Installed power optic system
- Push button control arm brake
- Adjustable handpiece drive air and mist control
- Selectable water source (tap or water bottle)
- Six bulb motion activated LED light with adjustable intensity and three axis light head movement (transformer included)
- 3+1 assistant instrumentation with integrated controls



## ADDITIONAL OPTIONS

- 4+1 handpiece configuration
- Hands free suction holder with extra HVE and tubing
- 8 bulb LED light and arm upgrade
- Flat screen display with mount (coming soon)

## A legend refined

### FDC50 PIVOT MOUNT PACKAGE

The 50 pivot mount package embodies maximum patient comfort and ergonomic workflow at an incredible value. You get the same versatile, reliable utility as our 38 Pivot Package but with our most advanced chair design and further improvements to the pivot mount delivery. With soft start/stop technology and updated chair movement, the hardest part of your job might be waking your patients up after you're done.

### CHAIR FEATURES

- Full hydraulic drive system with soft start/stop
- Advanced synchronized recline motion
- 50 degree locking chair swivel
- Locking swing down armrests, single handed operation
- Thin, extra wide aluminum backrest
- All metal frame and base construction
- Four position programmable chair memory
- Dual articulating, height adjustable headrest
- Plush Ultraleather upholstery




---

**FIRSTAR DENTAL CO.**

18915 E Valley Hwy  
Kent, WA 98032  
hello@firstardental.com  
[www.firstardental.com](http://www.firstardental.com)

---



Reliable. Purposefully designed. Incredible value. We got tired of choosing only two so we decided to fix that. For over 20 years, we've been designing and building dental equipment for every dentist, whether you're after the top of the line or something a bit more budget friendly that still checks all the boxes. We set out to make dental equipment easy because we know you want to focus on what you do best - dentistry. With Firstar, it's not wrong to love your dental equipment.

