



FDC50

SIDE MOUNT PACKAGE

DELIVERY STANDARD FEATURES

- Traditional delivery system with automatic control for three handpieces
- Pre-installed power optic system
- Push button control arm brake
- Selectable water source (tap or water bottle)
- Swiveling cuspidor with timed bowl rinse and automatic cup filler
- Integrated heated water system with quick disconnect outlet
- Six bulb motion activated LED light with adjustable intensity and three axis light head movement
- 3+1 assistant instrumentation with integrated controls



ADDITIONAL OPTIONS

- 4+1 handpiece configuration
- Hands free suction holder with extra HVE and tubing
- 8 bulb LED light and arm upgrade
- Flat screen display with mount



Our most powerful chair package

FDC50 SIDE MOUNT PACKAGE

Things that separate you from the ground - tires, mattresses, shoes, dental chairs - are always worth the premium. The FDC50 Side Mount Package pairs the refinement of our 50 chair with the versatility of our Side Mount Delivery. Driven by our latest soft/start hydraulic system, the 50 gives your patients the comfort of modern chair design while the Side Mount Delivery gives you all the tools and advanced water features to work efficiently and ergonomically.


CHAIR FEATURES

- Full hydraulic drive system with soft start/stop
- Advanced synchronized recline motion
- 50 degree locking chair swivel
- Locking swing down armrests, single handed operation
- Thin, extra wide aluminum backrest
- All metal frame and base construction
- Four position programmable chair memory
- Dual articulating, height adjustable headrest
- Plush Ultraleather upholstery



FIRSTAR DENTAL CO.

18915 E Valley Hwy
Kent, WA 98032
hello@firstardental.com
www.firstardental.com



Reliable. Purposefully designed. Incredible value. We got tired of choosing only two so we decided to fix that. For over 20 years, we've been designing and building dental equipment for every dentist, whether you're after the top of the line or something a bit more budget friendly that still checks all the boxes. We set out to make dental equipment easy because we know you want to focus on what you do best - dentistry. With Firstar, it's not wrong to love your dental equipment.

