

teknik dental®

Building since 1984

tegra MP200

tegra MP100S porcelain furnace was designed with years of experience for your quick, easy and superior quality production. It is our newest economical flagship superseding our history of proven classical furnace models.

Our modern design is built on top of the highest quality electronics, resistance and precision thermocouple to ensure a perfect firing experience...

Murat Parlar

Properties

Large colorful touchscreen
100 programs - 17 parameters
Perfect compatibility with all powders
Heaters enclosed in quartz spiral glass
Adjustable parameters during active cycle
Test Functions and detailed error messages



Parameters

Pre-Heat Temp. - °C
Pre-Heat Wait Time - min.sec
Lift Close Time - 10 Steps - min.sec
Dry Time - min.sec
Target Temp. - °C
Heat Rate - °C/min
Firing Time - min.sec
Vacuum Start - °C
Vacuum End - °C
Vacuum Wait Time - min.sec
Cooling Temp. (for Crome-Cobalt) - °C
Lift Open Time - min.sec

Settings

Vaccum Level - mbar Offset (Temp. Callibration) Stand By Temp. Language (Türkçe/English) Date - Time Test Functions

MP200_Brochure_EN 20230404